ISO 9001-2015,WHO, FDA, European CE Nortified & ISO-13485 Certified









ASE-2037 PERIMETER (Lister Type):

Automatically records field to 90 degree in all meridians. The Carrier mechanism is smooth, silent and invisible. The patient perceives only the test object to move. Spring Check, giving an audible click to each 15 degree of rotation of the arc in either direction Renders close scrutiny of dial plate unnecessary, and ensures accurate registration. The device can be thrown out of action if silent movement is desired.

Arc perimeter. Rotating crescent with sliding colour targets. Bulb holder with directional

shade. Adjustable double chin rest. Round chart holder (broken off). Integrally mounted case

of white, green and red colour targets.

The marker a needle is housed in a retractable tube, sprung to sheath the when not actually marking. A set of interchangeable and permanent washable test object is included with the instrument. The objects subtend angles of 4 degree, 2 degree, 1 degree, 1/2 degree and are in white, green and red blue objects in the last three sizes are also included. Supplies complete with Transformer to work on 220 Volt A/c mains With 100 charts.

Manufactured By: SUDARSHAN CHEMICALS India 324007 0744-2362345, 09414188018

CE, ISO 9001:2015, WHO, FDA & ISO-13485 CERTIFIED With first prize awarded at Haroti Udyog Mela year 2003





ASE-2036 PERIMETER (Smith Type): ALPHA

The apparatus comprise of a calibrates arc.

Revolving chart holder.

The object carrier which moves over the arc.

Contains five colors and five apertures of different diameter.

Calibrated scale.

At the back of the are arrangement is proved for fixing of chart which has concentric circles corresponding to the degree of arc.

Adjustable chin rest.

All fitted over a study base with recveptacle for keeping charts.

Table top use.

Complete with 500.

ASE-2037 PERIMETER (ListerType): ALPHA

Automatically records field to 90 degree in all meridians. The Carrier mechanism is smooth, silent and invisible. The patient perceives only the test object to move. Spring Check, giving an audible click to each 15 degree of rotation of the arc in either direction Renders close scrutiny of dial plate unnecessary, and ensures accurate registration. The device can be thrown out of action if silent movement is desired.

Arc perimeter. Rotating crescent with sliding color targets. Bulb holder with directional shade. Adjustable double chin rest. Round chart holder (broken off). Integrally mounted case of white, green and red color targets.

The marker a needle is housed in a retractable tube, sprung to sheath the when not actually marking. A set of interchangeable and permanent washable test object is included with the instrument.

The object subtend angles of 4 degree, 1 degree, 1/2 degree, 1/4 degree and are in white, green and red blue objects in the last three sizes are also included. Supplies complete with Transformer to work on 220 Volt A/C mains With 100 Charts.

Manufactured By: SUDARSHAN CHEMICALS India 324007 0744-2362345, 09414188018

CE, ISO 9001:2015, WHO, FDA & ISO-13485 CERTIFIED With first prize awarded at Haroti Udyog Mela year 2003